



Automated ESR System



Main features

- » Works directly off primary EDTA tubes
- » Consumables free
- » Westergren method
- » Sedimentation curve available
- » 32 results in only 25 minutes





Automated ESR System

FAST AND EASY



RELIABLE

- Analyzing method: AI machine learning algorithm
- Measuring principle: infrared photometry
- Integrated quality control
- Maintenance free
- Environmentally friendly: no waste fluids

CONNECTIVITY



Wiener Lab Switzerland S. A.
Av. Cardinal - Mermillod 36
1227 Geneva - Switzerland
marketing@wiener-lab.com

Facebook: Wiener lab.
Twitter: @Wiener_Lab
Instagram: @Wienerlabgroup
LinkedIn: Wiener lab Group



www.wiener-lab.com